


MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH A NEPHELOMETRIC COMMERCIAL PROCEDURE
- RANGE OF ANALYSIS.....0.0 – 500 NG/L
- GOOD PRECISION.....INTER-RUN: CVs < 1.6 %
.....INTRA-RUN: CVs < 3.3 %
- LIQUID REAGENTSREADY TO USE
- SELF-LIFE.....18 MONTHS
- PROZONE> 5500 NG/L
- LDD< 5.2 NG/L

GENERAL APPLICATION



PRECISION

(*) LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE FERRITIN CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 NG/ML STANDARD ONE

CONTROL	LOW	HIGH
MEAN NG/ML	182.7	287.2
SD NG/ML	2.9	4.5
CV %	1.6	1.6

WITHIN-RUN PRECISION (N=10)

CONTROL	LOW	HIGH
MEAN NG/ML	180.5	286.6
SD NG/ML	5.9	9.3
CV %	3.3	3.2

BETWEEN-RUN PRECISION (N=10)

SENSITIVITY*

CONTROL	LEVEL 1
N	20
LDD NG/ML	<5.2


SENSITIVITY

ORDER INFORMATION


ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST

- REF: K-9510T FERRITIN BUFFER 25 ML
- REF: S-9310T FERRITIN CALIBRATOR SET 6 x 1 ML
- REF: C-9310T FERRITIN CONTROL SET..... 2 x 1 ML


CALIBRATION CURVE



PROZONE



METHOD COMPARISON



LINEARITY

