


## MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH TWO COMMERCIAL PROCEDURES
- RANGE OF ANALYSIS .....0.0 – 80.0 MG/L
- GOOD PRECISION.....INTER-RUN: CVs < 3.67 %  
.....INTRA-RUN: CVs < 2.15 %
- LIQUID REAGENTS .....READY TO USE
- SELF-LIFE.....2 YEARS
- PROZONE .....> 160 MG/L
- LDD .....< 1.0 MG/L

## GENERAL APPLICATION



## PRECISION

CONTROL	LEVEL 1	LEVEL 2	LEVEL 3
MEAN MG/L	16.85	33.61	59.32
SD	0.36	0.28	0.88
CV %	2.15	0.82	1.49

WITHIN-RUN PRECISION (N=10)

CONTROL	LEVEL 1	LEVEL 2	LEVEL 3
MEAN MG/L	16.92	33.12	58.93
SD	0.46	0.69	2.17
CV %	2.70	2.09	3.67

BETWEEN-RUN PRECISION (N=10)

## SENSITIVITY\*

(\* ) LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE CYSTATIN C CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 MG/L STANDARD ONE

CONTROL	LEVEL 1
N	20
LDD MG/L	<1.0


SENSITIVITY

## ORDER INFORMATION


ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST

- OEM AND BULK PRESENTATIONS AVAILABLE



## CALIBRATION CURVE



## PROZONE



## METHOD COMPARISON



## LINEARITY

