


## MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH A NEPHELOMETRIC COMMERCIAL PROCEDURE
- RANGE OF ANALYSIS.....0.0 – 50.0 MG/L
- GOOD PRECISION.....INTER-RUN: CVs < 5.2 %  
.....INTRA-RUN: CVs < 5.5 %
- LIQUID REAGENTS .....READY TO USE
- SELF-LIFE .....2 YEARS
- PROZONE .....> 120 MG/L
- LDD .....< 0.18 MG/L

## GENERAL APPLICATION



## PRECISION

(\*) LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE C REACTIVE PROTEIN CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 MG/L STANDARD ONE

CONTROL	LEVEL 1	LEVEL 2	LEVEL 3
MEAN MG/L	1.80	5.57	41.16
SD MG/L	0.0943	0.1767	0.9252
CV %	5.20	3.17	2.25

WITHIN-RUN PRECISION (N=10)

CONTROL	LEVEL 1	LEVEL 2	LEVEL 3
MEAN MG/L	1.66	5.58	40.82
SD MG/L	0.0699	0.1135	2.2997
CV %	4.20	2.03	5.60

BETWEEN-RUN PRECISION (N=10)

## SENSITIVITY\*

CONTROL	LEVEL 1
N	20
LDD MG/L	<0.18


SENSITIVITY

## ORDER INFORMATION


ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST

- REF: K-9510T/HS ULTRASENSITIVE C REACTIVE PROTEIN BUFFER..... 30 ML  
ULTRASENSITIVE C REACTIVE PROTEIN LATEX..... 5 ML
- REF: S-9310T C REACTIVE PROTEIN CALIBRATOR ..... 1 ML
- REF: C-9310T/H C REACTIVE PROTEIN CONTROL HIGH ..... 1 ML
- REF: C-9310T/L C REACTIVE PROTEIN CONTROL HIGH ..... 1 ML


## CALIBRATION CURVE



## PROZONE



## METHOD COMPARISON



## LINEARITY

