

C REACTIVE PROTEIN (CRP)

REF: K-9510M

MONOREAGENT PROCEDURE

In vitro diagnostic reagents for the quantitative determination of C Reactive Protein (CRP) in serum by means of particle-enhanced turbidimetric immunoassay.

Diagnostic Relevance

C- reactive protein (CRP) is one of the acute phase proteins being synthesised by hepatocytes. The serum concentration of CRP increases during acute stages of diverse diseases associated with inflammation and tissue injury. Elevated CRP has been demonstrated in nearly all bacterial and fungal infections. In addition, it has been shown to be increased in other diseases as neoplasia, and rheumatic diseases as well as in major surgery.

The diagnosis usefulness of CRP is based on the velocity and on the magnitude of its increase. Serum concentrations are raised within hours of disease onset and the increase can be as much 2000-fold. A rapid fall of CRP levels indicates recovery.

Principle

This CRP test is based upon the reactions between C reactive protein (CRP) and latex-covalently bound antibodies against human CRP. CRP values are determined photometrically.

Reagents

Each CRP kit contains:

A.- Buffer - 45 mL of TRIS buffer, pH: 8.2, and 0.09 % sodium azide as preservative.

B.- Latex reagent - 5 mL of Polystyrene particles (0.5%) coated with goat antibodies anti-human-CRP serum in a glycine buffer (0.1 M, pH: 8.2), containing NaCl (0.15 M) and bovine serum albumin (0.5%). Preservative: Sodium azide 0.075%

C.- Calibrator - 1 mL. Human - based reference fluid. Preservative: sodium azide, 0.075 %. All raw materials of human origin used in the manufacture of this product showed no reactivity when tested for HBsAg, anti-HIV-1/2 and HCV with commercially available test methods. However, this product should be handled as though capable of transmitting infectious diseases

Reagent Preparation

Working Reagent is prepared with 1 part of Latex Reagent and 9 parts of Buffer Reagent. Prepare a fresh WR based on its workload. Shake gently the reagents before pipetting.

Calibration Curve and Controls

Analytical Range up to 90 mg/L.

Calibrator 1	100 µl of Biolatex CRP Calibrator*
Calibrator 2	100 µl of Calibrator 1 + 100 µl of Saline Solution
Calibrator 3	100 µl of Calibrator 2 + 100 µl of Saline Solution
Calibrator 4	100 µl of Calibrator 3 + 100 µl of Saline Solution
Calibrator 5	100 µl of Saline Solution

(* See values on the label or on the insert. Multiply by the appropriate factor.

For quality control use Biolatex Control or other suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Control must be assayed and evaluated as for patient samples.

Storage and Stability

Reagents in the original vial is stable to the expiration date on the vial label when capped and stored at +2 - +8°C. Immediately following the completion of an assay run, the reagent vial should be capped until next use in order to maximize curve stability. Once opened the reagent can be used within 1 month if stored tightly closed at +2 - +8°C after use. Do not freeze reagents.

The CRP latex reagent should have a white, turbid appearance free of granular particulate. Visible agglutination or precipitation may be a sign of deterioration, and the reagent should be discarded.

The CRP buffer reagent should be clear and colourless. Any turbidity may be sign of deterioration and reagent should be discarded.

WR is stable for up to one month at 4°C. **It is recommended that each Laboratory prepares a fresh Working Reagent based on its workload.**

Precautions

For in vitro diagnostic use only. Do not pipette by mouth. Reagents containing sodium azide must be handled with precaution. Sodium azide can

form explosive azides with lead and copper plumbing. Since absence of infectious agents cannot be proven, all specimens and reagents obtained from human blood should always be handled with precaution using established good laboratory practices. Disposal of all waste material should be in accordance with local guidelines. As with other diagnostic tests, results should be interpreted considering all other test results and the clinical situation of the patient.

Materials required

Spectrophotometric analyser. Saline solution. Controls.

Specimens

Fresh or deep frozen serum. CRP remain stable for 8 days at +2 to +8°C. If the test should be performed later, it is recommended to freeze the serum. Avoid successive freezing and thawing. Discard haemolysed or contaminated samples. Heavily lipaemic sera and turbid frozen serum samples must be cleared with a delipidating agent. Delipidation of samples do not affect the results of CRP in serum samples. The cleared patient serum sample must be used on the same day, as turbidity may reoccur.

Procedure

Wavelength	550nm		
Temperature	37°C		
Cuvette	1cm light path		
Measurement against distilled water blank.			
Bring the reagents at 37°C and pipette:			
	Calibrator	Sample	Blank
Calibrator	3 µl	---	---
Sample	---	3 µl	---
Distilled Water	---	---	3 µl
Work. Reagent	500 µl	500 µl	500 µl
Mix and measure absorbance immediately (A1) incubate 2 min (37°C), after incubation read absorbance (A2).			

Calculation

Plot the calibration curve and the sample concentration is obtained by interpolation the sample absorbance (A2-A1) in the calibration curve.

If is an one point calibration:

$$\frac{(A2-A1)_{\text{sample}} - (A2-A1)_{\text{blank}}}{(A2-A1)_{\text{calibrator}} - (A2-A1)_{\text{blank}}} \times \text{Calibrator Concentration}$$

Specific Performance Characteristics

As is well known, the analytical characteristics of a clinical chemistry reagent depend on both the reagents and the instrument used. Multicenter studies indicate important differences in analytical characteristics among similar instruments. Therefore, the data must be calculated by each instrument.

Reference Values

Values < 6 - 8 mg/L are within the normal range. Each laboratory should establish an expected range for the geographical area in which it is located.

Literature

- Hokama Y, Nakamura RM. C-Reactive protein: current status and future perspectives. J Clin Lab Anal 1987; 1: 15-27
- Hessian PA, Palmer DG. The presence and possible significance of C-Reactive protein in rheumatoid inflammation. J Rheumatol 1985 1985; 12:871-5.
- Okamura JM, Miyagi JM, Terada K, Hokama Y. Potential clinical applications of C-reactive protein. J Clin Lab Anal 1990; 4:231-5.
- Müller M, Mierau R, Wohltmann D. Interference of IgM rheumatoid factor with nephelometric C-reactive protein determinations. J Immunol Methods 1985; 80: 77-90.
- Young DS. Effects of Drugs on Clinical Laboratory Test. 5th Edition, AACC Press, 2000.
- Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two analytical methods. Application of linear regression procedures for method comparison studies. Part I. J Clin Chem Clin Biochem 1983; 21:709-20.
- Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 - 224

Fecha de Caducidad Expiry Date Date of Piromption Data di Scadenza Data Expiração Verwendbar bis	Temperatura de almacen Storage Temperature Temperatura de Conservazione Temperatura de Conservação Lagertemperatur	Número de Lote Lot Number Número de Lot Número di Lotto Número de Lote Chargen-Nr	Para Diagnóstico In Vitro For In Vitro Diagnostic Usage In Vitro Per Use Diagnostico In Vitro Utilizar em Diagnostico In Vitro In Vitro Diagnosticum	Número de catálogo Catalog Number Número de catalogue Numero di catalogo Número de catálogo Katalognummer	Conformidad Europea European Conformity Conformité aux normes européennes Conformità europea Conformidade com as normas europeias CE-Konformitätskennzeichnung	Fabricado por Manufactured by Fabriqué par Fabricado da Fabricado por Hergestellt	Reactivo Reagent Réactif Reagenti Reagenz	Calibrador Calibrator Calibratore Calibratore Kalibrator	Tampón Buffer Tampon Tampone Buffer Puffer	Liofilizado Lyophilised Lyophilisé Liofilo Liofilizado Lyophilisiert	Concentración Concentration Concentrazione Concentrazione Konzentration	Control Alto / Control Bajo Control High / Control Low Contrôle élevé / Contrôle Bas Controllo Alto / Controllo Basso Controllo Alto / Controllo Basso Kontrolle Hoch / Kontrolle Niedrig