

MICROALBUMIN (Albumin in Urine)

REF: K-9570M

MONOREAGENT PROCEDURE

In vitro diagnostic reagents for the quantitative determination of albumin in urine (MAU) by means of particle-enhanced turbidimetric immunoassay.

Diagnostic Relevance

Increased albumin excretion detectable only by sensitive immunoassay (microalbuminuria) has been used for some years as a predictor of incipient nephropathy and cardiovascular disease in diabetic patients. Microalbuminuria has also been associated with hypertension and increased risk of cardiovascular disease in non-diabetic patients. Microalbuminuria occurs in response to acute inflammatory conditions such as ischaemia, trauma, and thermal injury, surgery, pancreatitis, and inflammatory bowel disease. In many of these conditions albumin excretion increases within minutes or hours of the initiating stimulus and only lasts for 24 to 72 h. The degree of microalbuminuria is proportional to the severity of the inflammatory insult, is predictive of outcome, and is not associated with any other features of renal impairment.

Conventional dip-stick and acid precipitation tests for detecting protein in urine lack the sensitivity required to delineate this condition. Dip-stick may yield negative or trace results even when the albumin excretion rate is 10 or 20 times normal; and the rate must increase to 200 or 300 micrograms per minute ($\mu\text{g}/\text{min}$) before nephropathy becomes clinically apparent as persistent proteinuria. Interest in measuring subclinical elevations in the albumin excretion rate has focused on individuals with an already established diagnosis of diabetes or essential hypertension. Providing proper care is taken to minimise the influence of exercise and poor metabolic control of the albumin excretion rate, the urinary albumin level has proved to be an excellent predictor of the progression to overt nephropathy in both insulin-dependent and non-insulin dependent diabetes.

Principle

This MAU test is based upon the reactions between albumin and latex-covalently bound antibodies against human albumin. MAU values are determined photometrically.

Reagents

Each kit contains:

A - Buffer - 45 mL of phosphate buffer, pH: 8.5, < 0,1 % sodium azide as preservative.

B - Latex reagent - 5 mL of a suspension of latex microparticles covalently bound anti-albumin antibodies suspended in a neutral aqueous solution, and < 0,1 % sodium azide as preservative.

C - Dilution Buffer - 15 ml of buffer TRIS with 0.1% gelatine, pH: 7.0. Preservative : sodium azide < 0,1 %.

D - Calibrator - 1 mL. Human - based reference fluid. Preservative: sodium azide, 0.075 %. All raw materials of human origin used in the manufacture of this product showed no reactivity when tested for HBsAg, anti-HIV-1/2 and HCV with commercially available test methods. However, this product should be handled as though capable of transmitting infectious diseases

Reagent Preparation

Working Reagent is prepared with 1 part of Latex Reagent and 9 parts of Buffer Reagent. Prepare a fresh WR based on its workload. Shake gently the reagents before pipetting.

Calibration Curve and Controls

Analytical Range up to 250 mg/L.

Calibrator 1	100 μl of Biolatex MAU Calibrator*
Calibrator 2	100 μl of Calibrator 1 + 100 μl of Dilution Buffer
Calibrator 3	100 μl of Calibrator 2 + 100 μl of Dilution Buffer
Calibrator 4	100 μl of Calibrator 3 + 100 μl of Dilution Buffer
Calibrator 5	100 μl of Calibrator 4 + 100 μl of Dilution Buffer
Calibrator 6	100 μl of Dilution Buffer

(* See values on the label or on the insert. Multiply by the appropriate factor.

For quality control use Biolatex Control or other suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Control must be assayed and evaluated as for patient samples.

Storage and Stability

Reagents in the original vial is stable to the expiration date on the vial label when capped and

stored at +2 - +8°C. Immediately following the completion of an assay run, the reagent vial should be capped until next use in order to maximize curve stability. Once opened the reagent can be used within 1 month if stored tightly closed at +2 - +8°C after use. Do not freeze reagents. The MAU latex reagent should have a white, turbid appearance free of granular particulate. Visible agglutination or precipitation may be a sign of deterioration, and the reagent should be discarded. The MAU buffer reagent should be clear and colourless. Any turbidity may be sign of deterioration and reagent should be discarded. WR is stable for up to two weeks at 4°C. **It is recommended that each Laboratory prepares a fresh Working Reagent based on its workload.**

Precautions

For in vitro diagnostic use only. Do not pipette by mouth. Reagents containing sodium azide must be handled with precaution. Sodium azide can form explosive azides with lead and copper plumbing. Since absence of infectious agents cannot be proven, all specimens and reagents obtained from human blood should always be handled with precaution using established good laboratory practices. Disposal of all waste material should be in accordance with local guidelines. As with other diagnostic tests, results should be interpreted considering all other test results and the clinical situation of the patient.

Materials required

Spectrophotometric analyser. Controls.

Specimens

Use 12 or 24 hour collection. centrifuge urine specimens. Screen these specimens using an albumin test strip. If the result is negative (approx. below 300 mg/L), analyse the specimens undiluted. If the result is positive, dilute the specimen with specific protein sample diluent to obtain a concentration below 250 mg/L We recommend to dilute samples with dilution buffer.

Procedure

Wavelength	600 nm		
Temperature	37°C		
Cuvette	1 cm light path		
Measurement against distilled water blank.			
Bring the reagents at 37°C and pipette:			
	Calibrator	Sample	Blank
Calibrator	2 µl	---	---
Sample	---	2 µl	---
Distilled Water	---	---	2 µl
Work. Reagent	500 µl	500 µl	500 µl
Mix and measure absorbance immediately (A1) incubate 4 min (37°C), after incubation read absorbance (A2).			

Calculation

Plot the calibration curve and the sample concentration is obtained by interpolation the sample absorbance (A2-A1) in the calibration curve.

If is an one point calibration:

$$\frac{(A2-A1)_{\text{sample}} - (A2-A1)_{\text{blank}}}{(A2-A1)_{\text{calibrator}} - (A2-A1)_{\text{blank}}} \times \text{Calibrator Concentration}$$

Specific Performance Characteristics

As is well known, the analytical characteristics of a clinical chemistry reagent depend on both the reagents and the instrument used. Multicenter studies indicate important differences in analytical characteristics among similar instruments. Therefore, the data must be calculated by each instrument.

Reference Values

For timed overnight urine collections an albumin excretion rate greater than 20 µg/min is considered to abnormal. These data are to be interpreted as a guide. Each laboratory should establish its own reference intervals.

Literature

Winocour PH. Microalbuminuria Bmi 1992.1 304:1196-7 Marshall SM. Screening for microalbuminuria: which measurement, Diabetic Medicine 1991. 8: 706-11

Osherg Y et al. Effects of storage time and temperature on measurement of small concentrations of albumin in urine. Clin Chem 1990; 36:1428-30

Gosling P. Microalbuminuria: a sensitive indicator of non-renal disease?. Ann Clin Biochem 1995; 31439-41

Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two analytical methods. Application of linear regression procedures for method comparison studies. Part I. J Clin Chem Clin Biochem 1983; 21:709-20.

Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 - 224

Fecha de Caducidad Expiry Date Date de Péremption Data di Scadenza Data Expiração Verwendbar bis	Temperatura de almacen Storage Temperature Temperature de Conservation Temperatura de Conservazione Temperatura de Conservação Lagertemperatur	Número de Lote Lot Number Número de Lot Numero di Lotto Número de Lote Chargen-Nr	Para Diagnóstico In Vitro For In Vitro Diagnostic Usage In Vitro Per Uso Diagnostico In Vitro Utilizar em Diagnostico In Vitro In Vitro Diagnosticum	Número de catálogo Catalog Number Número de catalogue Numero di catalogo Número de catálogo Katalognummer	Conformidad Europea European Conformity Conformité aux normes européennes Conformità europea Conformidade com as normas europeias CE-Konformitätszeichen	Fabricado por Manufactured by Fabriqué par Fabricato da Fabbriçad por Herzestell	Reactivo Reagent Réactif Reagenti Reagent Reagens	Calibrador Calibrator Calibrateur Calibratore Calibrador Kalibrator	Tampón Buffer Tampon Tampone Buffer Puffer	Lyofilizado Lyophilised Lyophilisé Lyofiliato Lyophilizado Lyophilisiert	Concentración Concentration Concentration Concentrazione Concentração Konzentration	Control Alto / Control Bajo Control High / Control Low Contrôle élevé / Contrôle Bas Controllo Alto / Controllo Basso Controllo Alto / Controllo Basso Kontrolle Hooh / Kontrolle Niedrig