

## RHEUMATOID FACTOR (RF)

REF: K-9540M

### MONOREAGENT PROCEDURE

In vitro diagnostic reagents for the quantitative determination of Rheumatoid Factor (RF) in serum by means of particle-enhanced turbidimetric immunoassay.

### Diagnostic Relevance

The most consistent serological feature of rheumatoid arthritis is the increased concentration of autoantibodies directed against antigenic sites in the Fc region of human and animal IgG, namely rheumatoid factors (RFs) in the blood and joint fluid. The potential role of these factors in the pathogenesis of this disease has been studied extensively, with the finding that both environmental and genetic factors affect production of RF. RF determinations are clinically important for the diagnosis, prognosis, and assessment of therapeutic efficacy of rheumatoid arthritis. Although RFs may be found in all immunoglobulin classes, the RF most frequently detected in the laboratory is IgM type, present in about 75 - 80 % of adult patients with rheumatoid arthritis but in about 10 % of children with juvenile rheumatoid arthritis.

### Principle

This RF test is based upon the reactions between IgM-anti-IgG (RF) in patient's sample and latex-covalently bound human IgG. RF values are determined photometrically.

### Reagents

Each RF kit contains :

**A.- Buffer** - 45 mL of Phosphate buffer (0,05 M) pH: 7,0 containing NaCl (0,15 M), detergent and polyethyleneglycol.

Preservative : sodium azide < 1g/L

**B.- Latex reagent** - 7,5 mL of a suspension of latex microparticules covalently bound human IgG in a glycin buffer (0,1 M, pH: 8,2), containing NaCl (0,15 M) and bovine serum albumin (0,5%).

Preservative: Sodium azide 0,075%

**C.- Buffer Dil** - 15 ml of buffer TRIS, pH: 7.0.  
 Preservative : sodium azide < 1g/L

**D.- Calibrator** - 1 ml. Human - based reference fluid. Preservative: sodium azide, 0.075 %. All raw materials of human origin used in the manufacture of this product showed no reactivity when tested for HBsAg, anti-HIV-1/2 and HCV with commercially available test methods.

However, this product should be handled as though capable of transmitting infectious diseases

### Reagent Preparation

Working Reagent is prepared with 1 part of Latex Reagent and 6 parts of Buffer Reagent. Prepare a fresh WR based on its workload. (shake gently the reagents before pipetting).

**It is recommended that each Laboratory prepares a fresh Working Reagent based on its workload.**

### Calibration Curve and Controls

Analytical Range up to 140 IU/mL.

Calibrator 1	100 µl of Biolatex RF Calibrator*
Calibrator 2	100 µl of Calibrator 1 + 100 µl of Buffer Dil
Calibrator 3	100 µl of Calibrator 2 + 100 µl of Buffer Dil
Calibrator 4	100 µl of Calibrator 3 + 100 µl of Buffer Dil
Calibrator 5	100 µl of Buffer Dil

(\* See values on the label or on the insert. Multiply by the appropriate factor.

For quality control use Biolatex Control or other suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Control must be assayed and evaluated as for patient samples.

### Storage and Stability

Reagents in the original vial is stable to the expiration date on the vial label when capped and stored at +2 - +8°C. Immediately following the completion of an assay run, the reagent vial should be capped until next use in order to maximize curve stability. Once opened the reagent can be used within 1 month if stored tightly closed at +2 - +8°C after use. Do not freeze reagents.

The RF latex reagent should have a white, turbid appearance free of granular particulate. Visible agglutination or precipitation may be a sign of deterioration, and the reagent should be discarded.

The RF buffer reagent should be clear and colourless. Any turbidity may be sign of deterioration and reagent should be discarded.

## Precautions

For in vitro diagnostic use only. Do not pipette by mouth. Reagents containing sodium azide must be handled with precaution. Sodium azide can form explosive azides with lead and copper plumbing. Since absence of infectious agents cannot be proven, all specimens and reagents obtained from human blood should always be handled with precaution using established good laboratory practices. Disposal of all waste material should be in accordance with local guidelines. As with other diagnostic tests, results should be interpreted considering all other test results and the clinical situation of the patient.

## Materials required

Spectrophotometric analyser. Controls.

## Specimens

Serum specimens should be collected by venipuncture following good laboratory practices. RF remain stable for 72 hours at +2...+8°C. If the test should be performed later, it is recommended to freeze the serum. Heavily lipemic specimens, or turbid frozen specimens after thawing, must be clarified before the assay with a delipidating agent or by a high-speed centrifugation. Delipidation of samples do not affect the results of RF in serum samples. The cleared patient serum sample must be used on the same day, as turbidity may reoccur. Heat-inactivation of the sera is not necessary since C1q complement factor do not interfere in the assay.

## Procedure

Wavelength	600 nm		
Temperature	37°C		
Cuvette	1cm light path		
Measurement against distilled water blank.			
Bring the reagents at 37°C and pipette:			
	Calibrator	Sample	Blank
Calibrator	13 µl	---	---
Sample	---	13 µl	---
Distilled Water	---	---	13 µl
Work. Reagent	500 µl	500 µl	500 µl
Mix and measure absorbance immediately (A1) incubate 2 min (37°C), after incubation read absorbance (A2).			

## Calculation

Plot the calibration curve and the sample concentration is obtained by interpolation the sample absorbance (A2-A1) in the calibration curve.

If is an one point calibration:

$$\frac{(A2-A1)_{\text{sample}} - (A2-A1)_{\text{blank}}}{(A2-A1)_{\text{calibrator}} - (A2-A1)_{\text{blank}}} \times \text{Calibrator Concentration}$$

## Reference Values

Values <20 IU/ml are within the normal range.

This data has to be interpreted as a guide. Each laboratory should establish its own reference intervals.

## Specific Performance Characteristics

As is well known, the analytical characteristics of a clinical chemistry reagent depend on both the reagents and the instrument used. Multicenter studies indicate important differences in analytical characteristics among similar instruments. Therefore, the data must be calculated by each instrument.

## Literature

- Arnet FC, Edworthy SM, Bloch DA, McShane DJ, et al. The American Rheumatism Association 1987. Revised criteria for the classification of rheumatoid arthritis. *Arthritis Rheum* 1988; 31:315-24.
- Bartfield H. Distribution of rheumatoid factor activity in non rheumatoid states. *Ann NY Acad Sci* 1969; 168:30-40.
- Singer JM, Plotz CM. The latex fixation test. *Am J Med* 1956; 21:888-92.
- Moore TL, Dorner RN. Rheumatoid factors. *Clin Biochem* 1993; 26:75-84.
- Young DS. Effects of Drugs on Clinical Laboratory Test. 5th Edition, AACC Press, 2000.
- Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two analytical methods.
- Application of linear regression procedures for method comparison studies. Part I. *J Clin Chem Clin Biochem* 1983; 21:709-20.
- Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 - 224

Significados de los símbolos indicados en las etiquetas. Explanation of symbols used on labelling. Explication des symboles figurant sur les étiquettes. Spiegung der Symbole auf den Etiketten. Significado dos símbolos indicados nas etiquetas. Erläuterung der Symbole auf den etiketten.												
Fecha de Caducidad Expirate Date Date de Péremption Data di Scadenza Data Expiração Verwendbar bis	Temperatura de almacen Storage Temperature Temperatura de Conservation Temperatura de Conservazione Temperatura de Conservação Lagertemperatur	Número de Lote Lot Number Número de Lot Número di Lote Número de Lote Charge-Nr	Para Diagnóstico In Vitro For In Vitro Diagnostic Usage In Vitro Per Uso Diagnostico In Vitro Utilizar em Diagnostico In Vitro In Vitro Diagnosticum	Número de catálogo Catalog Number Número de catalogue Numero di catalogo Número de catálogo Katalognummer	Conformidad Europea European Conformity Conforme aux normes européennes Conformità europea Conformidade com as normas europeias CE-Konformitätskennzeichnung	Fabricado por Manufactured by Réalisé par Fabiicado de Fertiggestellt	Reactivo Reagent Réactif Reagent Reagente Reagenz	Calibrador Calibrator Calibrateur Calibratore Calibrador Kalibrator	Tampón Buffer Tampoon Tampone Buffer Puffer	Liofilizado Lyophilized Liofilisé Liofilizzato Liofilisiert	Concentración Concentration Concentrazione Concentração Konzentration	Control Alto / Control Bajo Control High / Control Low Controlle élevé / Contrôle Bas Controllo Alto / Controllo Basso Controllo Alto / Controllo Basso Kontrolle Hoch / Kontrolle Niedrig