

## $\alpha_1$ -MICROGLOBULIN ( $\alpha_1$ -m)

REF: L-9580T  
B-9480T

In vitro diagnostic reagents for the quantitative determination of alpha-1-microglobulin ( $\alpha_1$  - microglobulin) in urine by means of particle-enhanced turbidimetric immunoassay.

### Diagnostic Relevance

$\alpha_1$ - microglobulin ( $\alpha_1$ - m) is a low molecular weight glycoprotein (24,000- 33,000 g/mol) which was initially isolated from the urine of patients with renal tubular disorders in 1975. It is mainly synthesized in the liver and is widely distributed in various body fluids. Determination of  $\alpha_1$ - microglobulin in urine can be of aid in the diagnosis of tubular proteinuria. Detection of elevated concentrations of low molecular weight proteins in urine such as  $\alpha_1$ - microglobulin indicates tubular damage, which can occur in the course of advanced diabetic nephropathy, after heavy metal exposure or in the course of nephritis or others pathologies.

### Principle

In the development of alpha-1-microglobulin ( $\alpha_1$ -microglobulin) test, antibody to human  $\alpha_1$ - m bound to latex particles is brought into contact with  $\alpha_1$ - m in a sample. The increase in light scatter resulting from the agglutination reaction is proportional to the concentration of  $\alpha_1$ - m in the sample. The relationship between absorbance and concentration permits a multipoint calibration with a measuring range between 0 and 90 mg/L. Incubation takes place at 37°C. The assay can be performed on all instruments allowing turbidimetric measurements at 500 to 600 nm.

A range of standards are generated by successive dilution of the main standard in saline. The blank consists of saline.

### Reagents

A.- Buffer - TRIS buffer, pH: 7,2, containing detergents, polyethyleneglycol and <0.1 % sodium azide as preservative.

B.- Latex reagent - suspension of latex microparticules covalently bound antibodies to human  $\alpha_1$ - microglobulin suspended in a neutral aqueous solution, with < 0.1 % sodium azide as preservative.

### Precautions

For in vitro diagnostic use only. Do not pipette by mouth. Reagents containing sodium azide must be handled with precaution. Sodium azide can form explosive azides with lead and copper plumbing. Since absence of infectious agents cannot be proven, all specimens and reagents obtained from human blood should always be handled with precaution using established good laboratory practices. Disposal of all waste material should be in accordance with local guidelines. As with other diagnostic tests, results should be interpreted considering all other test results and the clinical situation of the patient.

### Materials required

Automatic analyzer. Saline solution. Calibrator. Controls.

### Storage and Stability

The  $\alpha_1$ - microglobulin reagents should be stored tightly capped at +2...+8°C when not in use. Do not freeze. Reagents in the original vials are stable to the expiration date on the vial label when capped and stored at +2...+8°C. Immediately following the completion of an assay run, the reagent vials should be capped until next use in order to maximize curve stability. Once opened the reagent can be used within 1 month if stored tightly closed at +2...+8°C after use. The  $\alpha_1$ - microglobulin buffer reagent should be clear and colourless. Any turbidity may be sign of deterioration and reagent should be discarded. The  $\alpha_1$ - microglobulin latex reagent should have a white, turbid appearance free of granular particulate. Visible agglutination or precipitation may be a sign of deterioration, and the reagent should be discarded.

### Specimens

Urine. All samples should be centrifuged prior to assay.  $\alpha_1$ - m remain stable in urine for 4 weeks at +2 to +8°C. If the test should be performed later, it is recommended to freeze the urine. Avoid successive freezing and thawing.

## Procedure

The reagents are ready to use as supplied. Latex reagent should be gently shaken (invert the recipient 3-4 times) before each use.

Volume R1/working reagent:	Volume R2/start reagent:	Volume sample:
250 µl	60 µl	2 µl
Step 1: mix R1 and R2, add sample and read 1st reading immediately after mixing.		
Step 2: 6 min after read 2nd reading.		
Wavelength: 500 nm		Incubation Time at 37° C: 7 min

\* Volume, time and wavelength are recommended. Adjust them depending of analyser features.

This reagent is intended to be used in clinical chemistry analysers. Adaptations for some of them are available.

## Calibration. Quality Control

Standardization: use Biolatex Calibrators. The method was standardized with reference to highly purified proteins preparation. The  $\alpha 1$ - m concentration of the Standard and Control is given on the label. Prepare the following dilutions of the standards using saline solution: 1; 1/2; 1/4; 1/8; 1/16, saline. The standard dilutions are to be used for measurement within 4 hours.

This curve is stored in memory by the analyser and recalled for later use. Calibration curves are stable for up to 14 days, after which a new curve must be generated. Additionally, recalibration must be performed whenever reagent lots are changed.

For quality control use BioLatex Control or other suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Control must be assayed and evaluated as for patient samples.

## Calculation

The turbidimetric analysers automatically calculate the  $\alpha 1$ - m concentration of each sample. Conversion: mg/l = µg/ml.

## Reference Values

Values < 10 mg/l (<14 mg/24h) are within the normal range. This data must be interpreted as a guide. Each laboratory should establish its own reference intervals.

## Automatic Analyzer

This product is performed for use it in turbidimetric automatic analysers or in manual procedures.

## Specific Performance Characteristics\*

As is well known, the analytical characteristics of a clinical chemistry reagent depend on both the reagents and the instrument used. Multicenter studies indicate important differences in analytical characteristics among similar instruments. Therefore, this data must be calculated by each instrument.

(\*) Analytical characteristics obtained in a single experiment in a Cobas-Mira plus analyser could be provided under demand.

## Literature

Berggard B, Ekstrom B and Akerstrom B.  $\alpha 1$ -microglobulin. Scand. J.Clin Lab. Invest. 1980;40; Suppl.154.

Kusano E, Suzuki M,Asano Y, Itoh Y, Takagi K and Kawai T. Human  $\alpha 1$ -microglobulin and its relationship to renal function. Nephron 1985;41;320-24

Straub, J.P; Baard M.A; du Jour H.A; Verplanke A.J.W and Herber F.M. The determination of  $\alpha 1$ -microglobulin by means of an automated latex immunoassay. Eur J Clin Chem clin Biochem 1995; 33;425-31.

Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 - 224.



BL, SL, C/Calahorra 4 - 6.26006 LO-Spain

May 2004

Significados de los símbolos indicados en las etiquetas. Explanation of symbols used on labelling. Explication des symboles figurant sur les étiquettes. Spiegung der Symbole auf den Etiketten. Significado dos símbolos indicados nas etiquetas. Erläuterung der Symbole auf den Etiketten.												
Fecha de Caducidad Expiry Date Date de Péremption Data di Scadenza Data Expiração Verwendbar bis	Temperatura de almacen Storage Temperature Temperature de Conservation Temperatura de Conservazione Temperatura de Conservação Lagertemperatur	Número de Lote Lot Number Número de Lot Numero di Lotto Número de Lote Chargen-Nr	Para Diagnóstico In Vitro For In Vitro Diagnostic Usage In Vitro Per Uso Diagnostico In Vitro Utilizar em Diagnostico In Vitro In Vitro Diagnosticum	Número de catálogo Catalog Number Número de catalogue Numero di catalogo Número de catálogo Katalognummer	Conformidad Europea European Conformity Conformité aux normes européennes Conformità europea Conformidade com as normas europeias CE-Konformitätskennzeichnung	Fabricado por Manufactured by Fabriqué par Fabricato da Fabricado por Hergestellt	Reactive Reagent Réactif Reagenti Reagente Reagens	Calibrador Calibrator Calibrateur Calibratore Calibrador Kalibrator	Tampón Buffer Tampon Tampone Buffer Puffer	Liofilizado Lyophilised Lyophilisé Liofilizzato Liofilizado Lyophilisiert	Concentración Concentration Concentration Concentrazione Concentração Konzentration	Control Alto / Control Bajo Control High / Control Low Contrôle élevé / Contrôle Bas Controlo Alto / Controlo Basso Controlo Alto / Controlo Baixo Kontrolle Hoch / Kontrolle Niedrig