


MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC MONOREAGENT METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH A COMMERCIAL PROCEDURE
- RANGE OF ANALYSIS.....0.0 – 250 MG/L
- GOOD PRECISION.....INTER-RUN: CVs < 3.7 %
.....INTRA-RUN: CVs < 4.5 %
- LIQUID REAGENTSREADY TO USE
- SELF-LIFE2 YEARS
- PROZONE> 500 MG/L
- LDD< 5.0 MG/L

GENERAL APPLICATION



(*) WR (WORKING REAGENT) IS PREPARED MIXING ONE PART OF REAGENT LATEX AND NINE PARTS OF REAGENT BUFFER

PRECISION

(*) LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE MICROALBUMIN CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 MG/L STANDARD ONE

LEVEL	1	2	3
N	10	10	10
CONCENTRATION MG/L	30	60	130
CV %	3.4	3.7	1.9

INTRA-ASSAY PRECISION

LEVEL	1	2	3
N	10	10	10
CONCENTRATION MG/L	30	60	160
CV %	4.5	2.4	4.5

INTER-ASSAY PRECISION

SENSITIVITY*

CONTROL	LEVEL 1
N	20
LDD MG/L	<5.00


SENSITIVITY

ORDER INFORMATION


ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST

- REF: K-9570M MICROALBUMIN BUFFER..... 45 ML
 MICROALBUMIN LATEX..... 5 ML
 MICROALBUMIN CALIBRATOR..... 1 ML
- REF: C-9370T MICROALBUMIN CONTROL..... 1 ML


CALIBRATION CURVE



PROZONE



METHOD COMPARISON



LINEARITY

