


MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH A NEPHELOMETRIC COMMERCIAL PROCEDURE
- RANGE OF ANALYSIS.....0.0 – 12.0 MG/L
- GOOD PRECISION.....INTER-RUN: CVs < 2.3 %
..... INTRA-RUN: CVs < 7.5 %
- LIQUID REAGENTSREADY TO USE
- SELF-LIFE.....2 YEARS
- PROZONE> 100 MG/L
- LDD< 0.2 MG/L

GENERAL APPLICATION



PRECISION

CONTROL	LEVEL 1	LEVEL 2	LEVEL 3
MEAN MG/L	1.8	5.8	11.5
SD MG/L	0.04	0.13	0.24
CV %	2.3	2.2	2.1

WITHIN-RUN PRECISION (N=10)

CONTROL	LEVEL 1	LEVEL 2	LEVEL 3
MEAN MG/L	1.8	5.8	11.5
SD MG/L	0.07	0.23	0.86
CV %	4.1	3.9	7.5

BETWEEN-RUN PRECISION (N=10)

SENSITIVITY*

(*) LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE β₂-MICROGLOBULIN CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 MG/L STANDARD ONE

CONTROL	LEVEL 1
N	20
LDD MG/L	<0.2

SENSITIVITY


ORDER INFORMATION

ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST


- REF: K-9530T β₂-MICROGLOBULIN BUFFER 30 ML
- REF: S-9330T β₂-MICROGLOBULIN LATEX 5 ML
- REF: C-9330T β₂-MICROGLOBULIN CALIBRATOR SET 1 ML
- REF: G-9330T/L β₂-MICROGLOBULIN CONTROL LOW..... 1 ML
- REF: G-9330T/H β₂-MICROGLOBULIN CONTROL HIGH 1 ML

biolatrix β_2 -MICROGLOBULIN


CALIBRATION CURVE



PROZONE



METHOD COMPARISON



LINEARITY

