

β₂-MICROGLOBULIN

REF: K-9530M

MONOREAGENT PROCEDURE

In vitro diagnostic reagents for the quantitative determination of β₂-microglobulin (β₂-m) in serum and urine by means of particle-enhanced turbidimetric immunoassay.

Diagnostic Relevance

β₂-microglobulin (β₂-m), an 11,800-Da protein consisting of 100 amino acid residues, first isolated from urine by Berggard and Bearn, is located on the surface of lymphocytes and other nucleated cells and was identified to be the light chain of the class I major histocompatibility complex.

The cell membrane turnover is the main source of β₂-m in serum, where its level rises when its production rate increases. Free β₂-m is filtered by the glomerulus and subsequently reabsorbed almost completely in the proximal tubular cells, where it is catabolized. Increased urinary excretion of β₂-m is a sensitive indicator of renal tubular disorders and has been used to detect early nephrotoxicity in patients treated with gentamicin and other nephrotoxic drugs.

Besides, renal insufficiency serum level of β₂-m was shown to be elevated in a variety of diseases including carcinomas and lymphoid tumors and inflammatory and autoimmune diseases such as Sjögren's syndrome and rheumatoid arthritis. Detection of elevated serum β₂-m has also been reported as a useful marker of acquired immune deficiency syndrome and in myeloma patients.

This test provides the advantage of being economical, rapid, precise, accurate and suitable for the analysis of large series of serum and urine specimens.

Principle

This β₂-m test is based upon the reactions between β₂-m in the sample and latex-covalently bound goat antihuman β₂-m antibodies. β₂-m values are determined photometrically.

Reagents

Each β₂-m kit contains:

A.- Buffer – 42.5 mL of TRIS buffer, pH: 7,2, containing detergents, polyethyleneglycol and 0,09 % sodium azide as preservative.

B.- Latex reagent – 8.5 mL of a suspension of latex microparticles covalently bound goat

antihuman β₂-m antibodies in a glycine buffer (0,1 M), containing NaCl (0,15 M) and bovine serum albumin (0,5%). Preservative: Sodium azide 0,075%.

C.- Calibrator – 1 mL. Human - based reference fluid. Preservative: sodium azide, 0.075 %. All raw materials of human origin used in the manufacture of this product showed no reactivity when tested for HBsAg, anti-HIV-1/2 and HCV with commercially available test methods. However, this product should be handled as though capable of transmitting infectious diseases

Reagent Preparation

Working Reagent is prepared with 1 part of Latex Reagent and 5 parts of Buffer Reagent. Prepare a fresh WR based on its workload. Shake gently the reagents before pipetting.

It is recommended that each Laboratory prepares a fresh Working Reagent based on its workload.

Calibration Curve and Controls

Analytical Range up to 20 mg/L.

Calibrator 1	100 µl of Biolatex β ₂ -m Calibrator*
Calibrator 2	100 µl of Calibrator 1 + 100 µl of Saline Solution
Calibrator 3	100 µl of Calibrator 2 + 100 µl of Saline Solution
Calibrator 4	100 µl of Calibrator 3 + 100 µl of Saline Solution
Calibrator 5	100 µl of Calibrator 4 + 100 µl of Saline Solution
Calibrator 6	100 µl of Saline Solution

(* See values on the label or on the insert. Multiply by the appropriate factor.

For quality control use Biolatex Control or other suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Control must be assayed and evaluated as for patient samples.

Storage and Stability

Reagents in the original vial is stable to the expiration date on the vial label when capped and stored at +2 - +8°C. Immediately following the completion of an assay run, the reagent vial should be capped until next use in order to maximize curve stability. Once opened the reagent can be used within 1 month if stored tightly closed at +2 - +8°C after use. Do not freeze reagents.

The β2-m latex reagent should have a white, turbid appearance free of granular particulate. Visible agglutination or precipitation may be a sign of deterioration, and the reagent should be discarded.

The β2-m buffer reagent should be clear and colourless. Any turbidity may be sign of deterioration and reagent should be discarded.

Precautions

For in vitro diagnostic use only. Do not pipette by mouth. Reagents containing sodium azide must be handled with precaution. Sodium azide can form explosive azides with lead and copper plumbing. Since absence of infectious agents cannot be proven, all specimens and reagents obtained from human blood should always be handled with precaution using established good laboratory practices. Disposal of all waste material should be in accordance with local guidelines. As with other diagnostic tests, results should be interpreted considering all other test results and the clinical situation of the patient.

Materials required

Spectrophotometric analyser. Saline solution. Controls.

Specimens

Serum specimens should be collected by venipuncture following good laboratory practices. Suitable assay specimens are human serum samples, as fresh as possible (stored up to 7 days at +2...+8°C) or deep-frozen. Any additional clotting or precipitation which occurs due to the freeze/thaw cycle should be removed by centrifugation prior to assay.

Lipemic specimens, or turbid frozen specimens after thawing, must be clarified before the assay by high-speed centrifugation (10 min at approx. 15.000 rpm).

Heat-inactivation of the serum samples can lead to a loss of the antigenic properties of the β2-m and must therefore be avoided.

Procedure

Wavelength	600 nm		
Temperature	37°C		
Cuvette	1cm light path		
Measurement against distilled water blank.			
Bring the reagents at 37°C and pipette:			
	Calibrator	Sample	Blank
Calibrator	3.5 µl	---	---
Sample	---	3.5 µl	---
Distilled Water	---	---	3.5 µl
Work. Reagent	500 µl	500 µl	500 µl
Mix and measure absorbance immediately (A1) incubate 5 min (37°C), after incubation read absorbance (A2).			

Calculation

Plot the calibration curve and the sample concentration is obtained by interpolation the sample absorbance (A2-A1) in the calibration curve.

If is an one point calibration:

$$\frac{(A2-A1)_{\text{sample}} - (A2-A1)_{\text{blank}}}{(A2-A1)_{\text{calibrator}} - (A2-A1)_{\text{blank}}} \times \text{Calibrator Concentration}$$

Reference Values

Values from 0,8 to 2,4 mg/L are within the normal range. Each laboratory should establish an expected range for the geographical area in which it is located.

Specific Performance Characteristics

As is well known, the analytical characteristics of a clinical chemistry reagent depend on both the reagents and the instrument used. Multicenter studies indicate important differences in analytical characteristics among similar instruments. Therefore, the data must be calculated by each instrument.

Literature

- F. Ortega y col.. Influencia pronóstica de la β2-microglobulina en el mieloma múltiple. Med. Clin. (Barc) 1992; 99: 645-648
- Lifson A.R, Hessel N A, Buchbinder S P, O'Malley P M, Barnhart L, Segal M, Katz M H, Holmberg S D. Serum β2-microglobulin and prediction of progression to AIDS in HIV infection. The Lancet 1992; 339: 1436-40.
- Weissel M, Scherak O, et al. Serum β2-microglobulin and SLE. Arthritis Rheum 1976; 19:986.
- Evrin P E, Paterson P A, Wide L, Berggard Y. Radioimmunoassay of β2-microglobulin in human biological fluids. Scand J clin Lab Invest 1971; 28: 439-43.
- Schardijn G, Stadius van Eps L. β2-microglobulin: Its significance in the evaluation of renal function. Kidney International 1987; 32: 635-41.
- Becker W. Die Standardisierung immunchemischer Plasmaprotein-Bestimmungen. Laboratoriumblätter 1980; 30: 25 - 32.
- Young DS. Effects of Drugs on Clinical Laboratory Test. 5th Edition, AACC Press, 2000.
- Sonderdruck aus DG Klinische Chemie Mitteilungen 1995; 26: 207 - 224

				REF			CAL	Buffer	LYOPH	Conc.	Control H / Control L	
Fecha de Caducidad Expiry Date Date de Péremption Data di Scadenza Data Expiração Verwendbar bis	Temperatura de almacen Storage Temperature Temperatura de Conservación Temperatura di Conservazione Temperatura de Conservação Lagertemperatur	Número de Lote Lot Number Número de Lot Número di Loto Número de Lote Chargen-Nr	Para Diagnóstico In Vitro For In Vitro Diagnostic Usage In Vitro Per Usu Diagnostico In Vitro Utilizar em Diagnostico In Vitro In Vitro Diagnostikum	Número de catálogo Catalog Number Número de catálogo Número di catalogo Número de catálogo Katalognummer	Conformidad Europea European Conformity Conformité aux normes européennes Conformità europea Conformidade com as normas europeias CE-Konformitätskennzeichnung	Fabricado por Manufactured by Réalisé par Fabbricatore Fabricado por Hergestellt	Reactivo Reagent Réactif Reagente Reagenz	Calibrador Buffer Calibrateur Calibratore Calibrador Kalibrator	Tampón Buffer Tampon Tampone Buffer	Liofilizado Lyophilized Lyophilisé Liofilizzato Liofilizzato Lyophilisiet	Concentración Concentration Concentrazione Concentração Konzentration	Control Alto / Control Bajo Control High / Control Low Controlle Elevé / Contrôle Bas Controllo Alto / Controllo Basso Controllo Alto / Controllo Basso Kontrolle Hoch / Kontrolle Niedrig