


MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC MONOREAGENT METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH A NEPHELOMETRIC COMMERCIAL PROCEDURE
- RANGE OF ANALYSIS0.0 – 500 NG/L
- GOOD PRECISION.....INTER-RUN: CVs < 3.38 %
.....INTRA-RUN: CVs < 3.88 %
- LIQUID REAGENTSREADY TO USE
- SELF-LIFE..... 18 MONTHS
- PROZONE> 7000 NG/L
- LDD< 5.0 NG/L

GENERAL APPLICATION



(* WR (WORKING REAGENT) IS PREPARED MIXING ONE PART OF REAGENT LATEX AND 2.5 PARTS OF REAGENT BUFFER

PRECISION

LEVEL	1	2	3
N	10	10	10
CONCENTRATION NG/ML	270	170	37
CV %	0.98	1.41	3.38

INTRA-ASSAY PRECISION

LEVEL	1	2	3
N	10	10	10
CONCENTRATION NG/ML	270	170	37
CV %	1.26	1.74	3.88

INTER-ASSAY PRECISION

SENSITIVITY*

(* LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE FERRITIN CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 NG/ML STANDARD ONE

CONTROL	LEVEL 1
N	20
LDD NG/ML	<5.00


SENSITIVITY

ORDER INFORMATION


ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST

- REF: K-956DM FERRITIN BUFFER37.5 ML
FERRITIN LATEX 15 ML
FERRITIN CALIBRATOR SET 1 ML
- REF: G-9360T FERRITIN CONTROL SET..... 2 x 1 ML


CALIBRATION CURVE




PROZONE



METHOD COMPARISON



LINEARITY




biolatrix **C** REACTIVE PROTEIN **MONOREAGENT**

MAIN FEATURES

- LATEX ENHANCED IMMUNOTURBIDIMETRIC MONOREAGENT METHOD
- EASY TO PROGRAM ON ANY TYPE OF AUTOMATED ANALYZERS
- GOOD CORRELATION WITH A NEPHELOMETRIC COMMERCIAL PROCEDURE
- RANGE OF ANALYSIS0.0 – 90.0 MG/L
- GOOD PRECISION.....INTER-RUN: CVs < 4.2 %
.....INTRA-RUN: CVs < 4.6 %
- LIQUID REAGENTSREADY TO USE
- SELF-LIFE.....2 YEARS
- PROZONE> 430 MG/L
- LDD< 2.0 MG/L

GENERAL APPLICATION



(*) WR (WORKING REAGENT) IS PREPARED MIXING ONE PART OF REAGENT LATEX AND NINE PARTS OF REAGENT BUFFER

PRECISION

CONTROL	LOW	MEDIUM	HIGH
CONCENTRATION MG/L	8.0	18.0	33.0
CV %	2.8	3.2	4.2

WITHIN-RUN PRECISION (N=10)

CONTROL	LOW	MEDIUM	HIGH
CONCENTRATION MG/L	8.0	18.0	33.0
CV %	3.9	4.1	4.6

BETWEEN-RUN PRECISION (N=10)

SENSITIVITY*

(*) LEAST DETECTABLE DOSE (SENSITIVITY) DETERMINED AS THE C REACTIVE PROTEIN CONCENTRATION NEEDED FOR A THREE DEVIATION RESPONSE OVER THE 0.0 MG/L STANDARD ONE

CONTROL	LEVEL 1
N	20
LDD MG/L	<2.00

SENSITIVITY

ORDER INFORMATION

ANOTHER PRESENTATIONS ARE AVAILABLE ON REQUEST

- REF: K-9510M C REACTIVE PROTEIN BUFFER 45 ML
C REACTIVE PROTEIN LATEX 5 ML
C REACTIVE PROTEIN CALIBRATOR 1 ML
- REF: C-9310T/H C REACTIVE PROTEIN CONTROL HIGH 1 ML
- REF: C-9310T/L C REACTIVE PROTEIN CONTROL HIGH 1 ML

INFO@BIOLATEX.COM

CALAHORRA, 4 – PAB. 6 – 26006 LOGROÑO (SPAIN)
TEL: (34) 941 27 00 47 FAX: (34) 941 27 10 77